(NL) CLAIRIVO ZWANGERSCHAPSTEST BASIS

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(EN) CLAIRIVO BASIC PREGNANCY TEST

INSTRUCTIONS FOR USE

Please read these instructions for use carefully before taking the test.

- For self-use
- Sample: urine
- Packaging: cassette

Application

Clairivo Basic Pregnancy Test is a self-use immunoassay test designed for rapid, visual and qualitative determination of human choriogonadotropin (HCG) in urine to aid in the early detection of pregnancy.

Summary

Human choriogonadotropin (HCG) is a glycoprotein hormone secreted by the developing placenta shortly after conception. During a normal pregnancy, HCG can be detected in the serum as early as 7 days after conception. The concentration of HCG continues to increase rapidly, often exceeding 100 mIU/ml at the time of the first missed menstrual period and reaching its maximum of 100,000 to 200,000 mIU/ml after a gestational age of 10-12 weeks. The detectability of HCG shortly after conception and the subsequent increase in concentration during the early development of the foetus makes this an excellent way of establishing early pregnancy.

Principle

Clairivo Basic Pregnancy Test is a chromatographic immunoassay (CIA) that enables rapid qualitative determination of HCG in urine. The membrane is coated in the test line area with anti- α HCG-binding antibodies and in the control line area with goat anti-mouse HCG. During testing, the urine sample may react with the coloured conjugate (mouse anti- β -HCG monoclonal antibody colloidal gold conjugate), which is pre-dried on the test strip. This mixture then moves chromatographically across the membrane by capillary action. A positive result is obtained if, at the site of the membrane test line, the specific, antibody-HCG-conjugate complex forms a pink line. The absence of this pink line at the site of the test line means that the result is negative. As the mixture continues to move across the membrane towards the immobilised goat anti-mouse HCG, a pink line will always form in the region of the control line, regardless of the presence or absence of HCG. The presence of this pink line indicates that 1) enough urine has been added, 2) the mixture has moved sufficiently and 3) the reagents have worked properly.

Composition

Clairivo Basic Pregnancy Test contains anti- β -HCG on colloidal gold particles and a combination of anti- α HCG on the membrane.

Warning and precautions

1. Read the instructions before performing this test.
2. Do not use after the indicated expiry date.
3. The test cassette should, until used, remain in the sealed container. Do not use if the packaging is damaged.
4. Do not use the test cassette more than once. Discard it after single use.
5. Do not touch the part of the membrane between the test lines.
6. Do not swallow the desiccant.

Storage

Store Clairivo Basic Pregnancy Test at between 2°C and 30°C in the sealed container, and not after the expiry date. Perform the test within 1 hour of opening the packaging. DO NOT FREEZE.


Collection of urine sample

1. Collect a fresh urine sample in a clean, dry container, which can be discarded after use.
2. Any urine sample can be used for a pregnancy test but a sample taken from the first urine of the day is best as the concentration of HCG in it is the highest.
3. A urine sample can be stored at room temperature for 8 hours. If the sample cannot be tested immediately, it may be stored at 2 - 8°C for 48 hours, or at -20°C for a long period. A urine sample should not be repeatedly frozen and thawed.
4. Bring the urine sample to room temperature before testing.

Test procedure

Allow the test cassette and urine sample to reach room temperature (15 - 30°C) before the taking the test. Do not open the packaging until the point of taking the test.

1. Remove the cassette from its sealed packaging and place it on a flat, clean and dry surface. Then use the cassette as soon as possible.
2. Use the pipette provided to aspirate a urine sample from the collection tray and put 3 drops into the collection circle.
3. Wait for the pink control line to become visible. Read the result within 5 minutes. Do not read the result after more than 5 minutes.

Interpreting the result

Pregnant: Two distinct pink lines are visible, one in the test zone (T) and one in the control zone (C).

NOTE: The intensity of the pink colour in the test zone (T) may vary depending on the concentration of HCG in the urine sample. Any degree of pink colouration in the test zone (T) should therefore be considered a positive result.

Not pregnant: Only one pink line is visible in the control zone (C). No clear pink line is visible in the test zone (T).

NOTE: No control line is visible.

NOTE: The most common causes of the control line not becoming visible are insufficient volume of urine sample and improper execution of the procedure. Check the procedure and repeat the test with a new test cassette. Contact your supplier if the problem persists.

Quality control

The test has its own control procedure. If a coloured line becomes visible in the control region (C), it indicates that the outcome of the internal control procedure is positive, that the test performed properly and the reagents worked well. If the background of the membrane remains visible, it indicates that the result of the internal control procedure is negative. When the test has been performed successfully and the reagents are working properly, the background will be clear and the result can be read properly.

Expected outcomes

Urine examination of healthy men and healthy, non-pregnant women will give a negative result. HCG will be present in the urine and serum of healthy, pregnant women. The amount of HCG will vary depending on the duration of pregnancy and per person. Clairivo Basic Pregnancy Test has an analytical sensitivity of 25 mIU/ml, and can demonstrate pregnancy as early as 1 day after missed menstruation.

Performance characteristics

Analytical sensitivity: Not less than 25 mIU/ml Analytical specificity: The result is negative for samples containing less than 500 mIU/ml hLH, 1000 mIU/ml hFSH and 1000 μ IU/ml hTSH.

Diagnostic sensitivity and diagnostic specificity:

This Clairivo Basic Pregnancy Test can detect a concentration of HCG when it is 25 mIU/ml or higher. 900 urine samples, known to be negative, were divided into 6 groups. HCG was added to each sample (150 per group) at a concentration of 0 mIU/ml, 5 mIU/ml, 10 mIU/ml, 25 mIU/ml, 50 mIU/ml or 5 IU/ml, calibrated according to the WHO 4th International Standard. Each group of samples was examined with Clairivo Basic Pregnancy Test. The results of this study were >99% consistent with the expected results.

Result	0 mIU/ml	5 mIU/ml	10 mIU/ml	25 mIU/ml	50 mIU/ml	5 IU/ml	Total
Positive	0	0	0	150	150	150	450
Negative	150	150	150	0	0	0	450
Total	150	150	150	150	150	150	900

Diagnostic sensitivity: 100% (450/450)

Diagnostic specificity: 100% (450/450)

Interference research: The following substances were added to urine samples without HCG and with HCG. None of the substances with the concentration examined affected the outcome of the determination.

• Acetaminophen:	0.2 mg/ml	Acetylsalicylic acid	0.2 mg/ml
• Ascorbic acid:	0.2 mg/ml	Atropine	0.2 mg/ml
• Caffeine:	0.2 mg/ml	Gentisic acid	0.2 mg/ml
• Glucose:	20 mg/ml	Haemoglobin	10 μ g/ml
• Tetracycline:	0.2 mg/ml		

Restrictions

1. Although it is not necessary to perform the test with the first morning urine, excessive fluid intake before the test should be avoided. The result may be "Not pregnant" if the urine sample is too diluted.
2. Package contents are for qualitative detection of HCG in urine only.
3. A sample with a low concentration of HCG may discolour over time. If a negative result is obtained while a pregnancy is suspected, a new sample should be collected after 48-72 hours and tested.
4. Fertility-promoting drugs that contain HCG may produce misleading results (these fertility-promoting drugs are usually administered by injection, and if the test is performed too soon after their administration, a pregnancy may be falsely established).
5. Other fertility-promoting treatments (such as using clomiphene citrate) and the use of painkillers and hormonal contraceptives (such as a contraceptive pill) will not affect the outcome.
6. HCG can be detected for several days to several weeks after childbirth, after a spontaneous abortion or after receiving injections with HCG.
7. Ectopic pregnancy, ovarian cysts, menopause and some very rare medical conditions can produce misleading results.
8. Although pregnancy is the most obvious reason for the presence of HCG in the serum and urine, an increase in HCG concentration without pregnancy has been reported in certain patients, such as in the case of trophoblastic disease and certain non-trophoblastic neoplasms.
9. Alcohol may affect the result of the test. It is recommended not to perform the test after drinking alcohol.
10. This test gives the presumptive diagnosis of "pregnancy". A pregnancy should only be diagnosed by a doctor after an assessment of all clinical and laboratory determinations.
11. The user with eye diseases such as color blindness or color weakness may interpret the result wrongly.

Questions & Answers

1. Q: How soon after first suspecting that I am pregnant can I take the test?
Answer: You can take the test as early as the first day after your missed period.
2. Q: Should I take the test using the first urine of the day?
Answer: You can take the test at any time of the day. The first morning urine is usually the most concentrated urine of the day and therefore contains the most HCG.
3. Q: How do I know I performed the test correctly? Answer: A coloured line becoming visible in the control area (C) indicates that you performed the test correctly and that enough urine was absorbed by the strip.
4. Q: Is the test considered valid after more than 5 minutes?
Answer: No. The result of the test must be read within 5 minutes. A positive result should not change within a few days, but a negative result can change to a false-positive result within minutes of the end of the test period, which is then not a valid result.

Bibliography


- [1] Robert D. Nerenz, Haowei Song, and Ann M. Gronowski.Screening Method to Evaluate Point-of-Care Human Chorionic Gonadotropin (hCG) Devices for Susceptibility to the Hook Effect by hCG β -Core Fragment: Evaluation of 11 Devices. Clinical Chemistry 60:4:667-674 (2014)
- [2] Berger P, Paus E, Hemken PM, Sturgeon C, Stewart WW, Skinner JP, Harwick LC, Saldana SC, Ramsay CS, Rupprecht KR, Olsen KH, Bidart JM,Stenman UH. Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-Workshop.Tumor Biol. (2013) 34:4033-4057

	Do not re-use		Use-by date		Batch code
	Manufacturer		Store at 2-30°C		CE mark
	In vitro diagnostic medical device		Consult instructions for use		Date of manufacture
	Contains sufficient for <n>- tests		Authorized representative in the European Union		Importer

	Wellkang Ltd. Enterprise Hub, NW Business Complex, 1 Beraghmore Rd. Derry, BT48 8SE, N. Ireland
--	--

	Core Technology Co., Ltd. Room 100, C Building, No. 29 Life Park Rd, Changping District, Beijing 102206, China
--	--

	Dayes BV. Exa 32, 6902 KH Zevenaar The Netherlands
--	--

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Littérature
 [1] Robert D. Nerenz, Haowei Song, and Ann M. Gronowski.Screening Method to Evaluate Point-of-Care Human Chorionic Gonadotropin (hCG) Devices for Susceptibility to the Hook Effect by hCG βCore Fragment: Evaluation of 11 Devices. Clinical Chemistry 60:4667–674 (2014)
 [2] Berger P, Paus E, Hemken PM, Sturgeon C, Stewart WW, Skinner JP, Harwick LC, Saldana SC, Ramsay CS, Rupprecht KR, Olsen KH, Bidart JM,Stenman UH. Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7workshop.Tumor Biol. (2013) 34:4033–4057

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Literatur
 [1] Robert D. Nerenz, Haowei Song, and Ann M. Gronowski.Screening Method to Evaluate Point-of-Care Human Chorionic Gonadotropin (hCG) Devices for Susceptibility to the Hook Effect by hCG βCore Fragment: Evaluation of 11 Devices. Clinical Chemistry 60:4667–674 (2014)
 [2] Berger P, Paus E, Hemken PM, Sturgeon C, Stewart WW, Skinner JP, Harwick LC, Saldana SC, Ramsay CS, Rupprecht KR, Olsen KH, Bidart JM,Stenman UH. Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7workshop.Tumor Biol. (2013) 34:4033–4057

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